



Display Compatibility: 26" Wide ( Weight)  
 Supported Weight: 13-18 Kg  
 VESA standard: 75 / 100 mm

Up to 26" Wide ( Weight)  
 75/100 mm  
 16-22 Kg

Do not hesitate to contact us if you have any questions about this product.

**Warranty: 5 years**

000.991.114-5

Graphite Grey and White

On Measure

= Depending on the elements

991.114-5

**PH Medical Monitor Arm (High Position) height adjustable 300 + 240 mm, wall mounting on vertical slide rail**

This PH (High Position) Horizontal Medical Arm 300 + 240 mm with gas spring height adjustment allows the LCD screen to be perfectly fixed, oriented and adjusted. It is coated with an anti-microbial agent offering hygiene and cleaning down to the last detail. Its modern design has been specially designed for healthcare environments.

Technical specifications:

- \* Integrated cable channel
  - \* Equipotential bonding
  - \* Optional Vertical Wall Slide Rail Mounting, including several lengths are available (480 mm - 720 mm - 960 mm - 1200 mm). These wall rails can be easily cut to size you want. All cables also pass through the wall rail. (See Option below)
  - \* This PH Horizontal 300 +240 mm gas spring arm simply slips into the wall rail and then adjusts to the desired height and fixed. It has a safety height adjustment stop button, which allows the arm to remain in the chosen position.
  - \* Installation specific to the high position (PH)
  - \* The head of this arm complies with the screen mounting standard: VESA MIS-D 75 mm and 100 mm.
  - \* Compatible with screens whose weight does not exceed 16-22 Kg.
  - \* Touch screen stability. \* Screen tilt: 50° down and 100° up
  - \* Screen rotation: 110° right, 110° left.
  - \* Arm rotation: 105° to the right, 105° to the left.
  - \* Height adjustment: 45° down, 45° up \* Total length 806 mm
  - \* Compatible with multimedia terminals (TMM)
  - \* Colour RAL 7024 (Graphite grey) and White \* Warranty: 5 years
  - \* All our medical arms are in conformity: CE, ROHS, Medical Grade, Regulations MDD 93/42 ECC
- Note: This arm can also be installed on one of our medical columns or trolleys.