

Part Number: 991.540-1—Presentation

300 + 240 mm height adjustable medical monitor arm, column


mounting












Display Compatibility: 26" Wide (Weight) Supported Weight: 3-8 Kg VESA standard: 75 / 100 mm





Do not hesitate to contact us if you have any questions about this product.



Warranty: 5 years



000.991.540-1



Graphite Grey and White



On Measure



= Depending on the elements



991.540-1

## 300 + 240 mm height adjustable medical monitor arm, column mounting

This 300 + 240 mm Horizontal Medical Arm with gas spring height adjustment allows the LCD screen to be perfectly fixed, oriented and adjusted. It is coated with an anti-microbial agent offering hygiene and cleaning down to the last detail. Its modern design has been specially designed for healthcare environments

## **Technical specifications:**

- \* Integrated cable channel
- \* Equipotential bonding
- \* Optional mounting on worktop or floor or ceiling column, several Heights are available:

(300 mm - 500 mm - 1500 mm - 2000 mm)

These columns can be easily cut to the desired size. All cables also pass through the Column. (See Option below)

- \* This 300 + 240 mm Horizontal arm with gas spring simply slips in in the Column and then adjust to the desired height and fixed. It has a safety stop button for the height adjustment, which allows the arm to remain in the chosen position.
- \* The head of this arm complies with the screen mounting standard: VESA MIS-D 75 mm and 100 mm.
- \* Compatible with screens whose weight does not exceed 3-8 Kg.
- \* Touch screen stability. \* Screen tilt: 50° down and 100° up
- \* Screen rotation: 110° right, 110° left.
- \* Arm rotation: 105° to the right, 105° to the left.
- \* Height adjustment: 45° down, 45° up \* Total length 806 mm
- \* Compatible with multimedia terminals (TMM)
- \* Colour RAL 7024 (Graphite grey) and White \* Warranty: 5 years
- \* All our medical arms are in conformity:

CE, ROHS, Medical Grade, Regulations MDD 93/42 ECC.

Note: This arm can also be installed on one of our medical columns or trolleys.